

# BIO II

## Biosafety Cabinet Type II



**Triple Protection :**  
Operator  
Product  
Environment

*Complies with standard:*  
**EN 12469 & NSF 49**

- ✓ compact size
- ✓ Automatic flow regulation ECM
- ✓ flip window
- ✓ UV Lamp

★ **FLAGSHIP PRODUCT**

The BIO II is adapted to handling of agents that require biosafety level class 1, 2 or 3.

	BIO II 9		BIO II 12	
	internal	Overall	internal	Overall
Width	919	922	1226	1229
Height	internal : 500		Overall : 1080	
Depth	internal: 495		Overall : 640	
Weight (Kg)	120		150	
Power (W)	1500			

