

# PPS

## Biosafety Cabinet Type II



**Triple protection :**  
Operator  
Product  
Environment

*Complies with standard:*  
**NF EN 12469**

- ✔ Electric sash window
- ✔ Bag in Bag out System
- ✔ Maintenance from the front
- ✔ 3 HEPA filters

★ **FLAGSHIP PRODUCT**

Microbiological safety stations type II are designed to protect the operator, product and the environment from any airborne contamination.

The PPS is suitable for handling pathogenic germs of class 1, 2 and 3.

The PPS can be fitted with a marble for accurate weighing at (0,01 mg)



	PPS 9	PPS 12	PPS 16	PPS 18
Useful Width	921	1226	1531	1836
Total width	971	1276	1581	1886
Height	internal : 650		overall : 2200	
Depth	Useful: 616		Overall: 786	
Weight (kg)	270	330	420	490
Power (KW)	1,4	1,4	1,5	1,5