

UNIVERS MV

Vertical flow cabinet - Class ISO 5



Protection :
Product

Complies with standard:
ISO 14644

- ✓ Manual regulation
- ✓ Economical solution
- ✓ Glass volume
- ✓ Flip window

The cabinet ISO class 5, are designed to protect the product thanks to a vertical flow of sterile air.

The enclosure is easy to install and can be installed by you.

| | UNIVERS MV 9 | |
|-------------|--------------|---------|
| | internal | Overall |
| Width | 905 | 925 |
| Height | 700 | 1030 |
| Depth | 592 | 600 |
| Weight (kg) | 30 | |