

Accessories

ATP Thermal Printer



Technical Specifications

Model	ATP Thermal Printer
Power Supply	External power adapter 110/220 VAC - 8.5 VDC OUT
Baud Rates	Selectable from 2400 to 115200
Data Bits	8
Parity	None
Operating Temperature	0°C to + 45°C
Housing	ABS plastic housing
Overall Dimensions (WxDxH)	4.5" x 7.2" x 3.5" / 114mm x 184mm x 90mm
Net Weight	1.7lb / 0.76kg with paper roll installed

Features

- Interface capabilities: serial port
- Support for up to 26 international character sets
- Print size 24 x 24 dot
- Column size 384 dots per line
- Print speed 50mm / sec
- 57.5mm paper width
- Selectable baud rates
- LED indicators to show power, error and when out of paper

Accessories

- No. 3126011263 Thermal paper
- No. 3126011281 Thermal paper (10pk)



PMB Temperature Calibration kit



Cushioned Carry Case

Technical Specifications

Model	PMB Temperature Calibration Kit
Power Supply	Battery - 1.5v DC
Display Value	0.1°C
Operating Temperature	-50°C to 300°C

PMB Moisture Sample Pans



Available in packs of 250

Technical Specifications

Model	PMB Sample Pans
Pan Size	100mm ø Top - 100mm - wide Depth - 8mm - height Bottom - 93mm - wide
Pan Material	Aluminium

