



CAT. 90008700

# DOT

All that counts  
Colony Counter



labopolis 

[www.labopolis.com](http://www.labopolis.com)

## Ordering Information

CAT. 90008700 DOT Colony Counter with 120mm magnifying glass with light.

CAT. 90008701 DOT Colony Counter with 100mm magnifying glass without light.

## Accessories and Consumables

1. Clickable permanent marker  
CAT. 900011459


2. Mixer Grid 1/6 cm<sup>2</sup>  
CAT. 9000613

3. 150mm centering device  
CAT. 9000615

4. 60mm centering device  
CAT. 9000606

5. Spiral Grid  
CAT. 9000625

6. Stainless contact pointer  
CAT. 900011460



## Technical data

Counter: 4 digits, 9 mm LED-Display. Maximum count 9999

Inputs: Clickable permanent marker/ Stainless contact pointer (optional).

Outputs: USB Human Interface Device (HID) emulating standard PC keyboard output receivable by PC running keyboard compatible applications (Microsoft Windows®).

USB Virtual COM port (VCP) reserved for IUL, S.A. Technical Service purposes.

Upper light power (12V DC, 500mA)

Magnifying glasses: Biconvex, 100/120 mm Ø, 4 Diopters, DIN x 2

Lightning: Bottom light, 4 selectable intensities, On/Off switching powered from the instrument.

Upper light (depending of the glass magnifier ordered), On/Off switching powered from the instrument.

Specifications subjects to change without notice.

Mains: 100/240 V, 50/60 Hz

Power: 50W

Dimensions (WxDxH):

211x320x335 mm / 8.3x12.6x13.2 in (Cat. No. 90008700)

211x320x395 mm / 8.3x12.6x15.5 in (Cat. No. 90008701)

Weight:

3,2 Kg / 7,0 lb (Cat. No. 90008700)

2,7 Kg / 5,95 lb (Cat. No. 90008701)

For more information please visit us at:

[www.STabaTe.com](http://www.STabaTe.com)

>34ABA>Ei Ež>ž>a` VdW# " Ž@ShW) - \$\*\*"' 3 1S'VW ž? SvqV (Spain)

[ Xa@STabaTe.com · T. +34 9# \$( ' + ' ( "

