

OPTIGEL

Horizontal flow cabinet - Class ISO 5



Protection :
Product

Complies with standard:
ISO 14644

- ✔ Opens completely
- ✔ ECM regulation
- ✔ Plexiglas side windows
- ✔ Stainless steel work surface

The enclosure Type OPTIGEL consist of a blowing and filtration unit. Optimum handling protection thanks to a sterile air flow of class ISO 5 .

	OPTIGEL 9		OPTIGEL 12		OPTIGEL 18	
	internal	Overall	internal	Overall	internal	Overall
Width	920	925	1225	1250	1820	1840
Height	INTERNAL : 550 OVERALL : 813					
Depth	INTERNAL : 550 OVERALL : 911					
Height base	1740					
Weight (kg)	80		100		150	